

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

### 1 Identification of the substance/mixture and of the company/undertaking

- **Product identifier**
- **Trade name: EpoFix Hardener 130ml**
- **Article number:** AGB8790H
- **CAS Number:**  
112-24-3
- **EC number:**  
203-950-6
- **Index number:**  
612-059-00-5
- **Relevant identified uses of the substance or mixture and uses advised against**  
For embedding and impregnation of materialographic specimens.
- **Application of the substance / the preparation:** No further data available.
- **Details of the supplier of the safety data sheet**
- **Supplier.**  
Agar Scientific  
Elektron Technology UK Ltd  
Parsonage Lane  
Stansted CM24 8GF  
United Kingdom  
sales@agarscientific.com  
Tel: +44 (0) 1279 813 519
- **Further information obtainable from:** Technical Support
- **Emergency telephone number:** 24 hours: +44 (0)1856 407333

### 2 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**



GHS05 corrosion

Skin Corr. 1B      H314 Causes severe skin burns and eye damage.



GHS07

Acute Tox. 4      H312 Harmful in contact with skin.

Skin Sens. 1      H317 May cause an allergic skin reaction.

Aquatic Chronic 3      H412 Harmful to aquatic life with long lasting effects.

- **Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**  
The substance is classified and labelled according to the CLP regulation.
- **Hazard pictograms**



GHS05    GHS07

- **Signal word** Danger
- **Hazard-determining components of labelling:**  
triethylenetetramine

(Contd. on page 2)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

**Trade name: EpoFix Hardener 130ml***(Contd. of page 1)***· Hazard statements**

- H312 Harmful in contact with skin.
- H314 Causes severe skin burns and eye damage.
- H317 May cause an allergic skin reaction.
- H412 Harmful to aquatic life with long lasting effects.

**· Precautionary statements**

- P101 If medical advice is needed, have product container or label at hand.
- P102 Keep out of reach of children.
- P103 Read label before use.
- P260 Do not breathe dust/fume/gas/mist/vapours/spray.
- P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 Immediately call a POISON CENTER or doctor/physician.
- P405 Store locked up.
- P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

**· Other hazards**

Corrosive. Prolonged contact causes serious eye and tissue damage. Be aware that symptoms of lung oedema (shortness of breath) may develop up to 24 hours after exposure. The product contains a sensitising substance which may provoke an allergic reaction among sensitive individuals.

**· Results of PBT and vPvB assessment**

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

### 3 Composition/information on ingredients

- **Chemical characterization: Substances**
- **CAS No. Description**  
112-24-3 triethylenetetramine
- **Identification number(s)**
- **EC number:** 203-950-6
- **Index number:** 612-059-00-5

### 4 First aid measures

**· Description of first aid measures****· General information:**

Immediately remove any clothing soiled by the product.  
Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

**· After inhalation:**

Move into fresh air and keep at rest. Be aware that symptoms of lung oedema (shortness of breath) may develop up to 24 hours after exposure. Immediately call an ambulance.

**· After skin contact:**

Remove contaminated clothes and rinse skin thoroughly with water. In case of eczema or other skin disorders: Seek medical attention and bring these instructions.

**· After eye contact:**

Immediately flush with plenty of water. Remove any contact lenses and open eyes wide apart. Call an ambulance and continue flushing during transportation to hospital. Bring these instructions.

**· After swallowing:**

Immediately rinse mouth and drink plenty of water. Keep person under observation. If person becomes uncomfortable seek hospital and bring these instructions. Do not induce vomiting.

*(Contd. on page 3)*

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

**Trade name: EpoFix Hardener 130ml***(Contd. of page 2)*

- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed**  
See section 11 for more detailed information on health effects and symptoms.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
Use fire-extinguishing media appropriate for surrounding materials.
- **Special hazards arising from the substance or mixture**  
When heated and in case of fire, very toxic nitrogen oxides may be formed.
- **Advice for firefighters**
- **Protective equipment:**  
Selection of respiratory protection for fire fighting: follow the general fire precautions indicated in the workplace.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**  
Avoid contact with skin and eyes.
- **Environmental precautions:** Do not discharge into drains, water courses or onto the ground.
- **Methods and material for containment and cleaning up:**  
Absorb spillage with non-combustible, absorbent material.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Avoid contact with skin and eyes. Immediately change drenched clothing. Observe good chemical hygiene practices.  
Local exhaust is recommended. Provide easy access to water supply and eye wash facilities.  
Work practice should minimise contact.
- **Information about fire - and explosion protection:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**  
Store in a cool and well-ventilated place.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep container tightly sealed.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.
- **Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:** Not required.

*(Contd. on page 4)*

# Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

**Trade name: EpoFix Hardener 130ml**

(Contd. of page 3)

- **Additional information:** No occupational exposure limit assigned.
- **Exposure controls**  
Mix and prepare in a place with efficient exhaust ventilation. Local exhaust is recommended. An eye wash bottle must be available at the work site.
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Keep away from foodstuffs, beverages and feed.  
Remove all soiled and contaminated clothing  
Wash hands before breaks and at the end of work.  
Avoid contact with the eyes and skin.
- **Respiratory protection:**  
In case of inadequate ventilation and work of brief duration, use suitable respiratory equipment.  
Use respiratory equipment with gas filter, type K.
- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.  
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

- **Material of gloves**  
Use disposable gloves protecting against epoxy along with cotton gloves closest to the skin. Laminate (PE/EVOH) gloves are recommended.
- **Penetration time of glove material**  
The most suitable glove must be chosen in consultation with the gloves supplier, who can inform about the breakthrough time of the glove material.
- **Eye protection:**



Tightly sealed goggles

Face shield.

- **Body protection:**  
Wear suitable protective clothing as protection against splashing or contamination.
- **Risk management measures**  
Personal protection equipment should be chosen according to the CEN standards and in discussion with the supplier of the personal protective equipment.

## 9 Physical and chemical properties

### · Information on basic physical and chemical properties

#### · General Information

#### · Appearance:

Form:	Liquid
Colour:	Yellow
Odour:	Amine-like
Odour threshold:	Not determined.

· pH-value: 10.7 (ISO 8975)

#### · Change in condition

    Melting point/Melting range: -35 °C

(Contd. on page 5)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

**Trade name: EpoFix Hardener 130ml**

(Contd. of page 4)

<b>Boiling point/Boiling range:</b>	270-300 °C
· <b>Flash point:</b>	129 °C
· <b>Flammability (solid, gaseous):</b>	Not applicable.
· <b>Ignition temperature:</b>	335 °C
· <b>Decomposition temperature:</b>	Not determined.
· <b>Self-igniting:</b>	Not determined.
· <b>Danger of explosion:</b>	Product does not present an explosion hazard.
· <b>Explosion limits:</b>	
<b>Lower:</b>	Not determined.
<b>Upper:</b>	Not determined.
· <b>Vapour pressure at 20 °C:</b>	<0.1 hPa
· <b>Density at 20 °C:</b>	0.982 g/cm <sup>3</sup>
· <b>Relative density</b>	Not determined.
· <b>Vapour density</b>	Not determined.
· <b>Evaporation rate</b>	Not determined.
· <b>Solubility in / Miscibility with water:</b>	Fully miscible.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
<b>Dynamic at 25 °C:</b>	20-30 mPas
<b>Kinematic:</b>	Not determined.
· <b>Other information</b>	No further relevant information available.

## 10 Stability and reactivity

- **Reactivity**
- **Chemical stability** Stable under normal temperature conditions.
- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No data available.
- **Conditions to avoid** Heating. Direct sunlight.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:**  
When heated and in case of fire, very toxic nitrogen oxides are formed.

## 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

- **LD/LC50 values relevant for classification:**

### 112-24-3 triethylenetetramine

Oral	LD50	2500 mg/kg (rat)
------	------	------------------

Dermal	LD50	805 mg/kg (rabbit)
--------	------	--------------------

- **Primary irritant effect:**

- **on the skin:**

Corrosive. Prolonged contact causes serious tissue damage. Risk of sensitisation or allergic reactions among sensitive individuals.

(Contd. on page 6)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

**Trade name: EpoFix Hardener 130ml**

*(Contd. of page 5)*

- **on the eye:** Corrosive. Immediate first aid is necessary.
- **Sensitization:** Sensitization possible through skin contact.
- **Other information (about experimental toxicology):**  
 Inhalation: Not relevant at normal room temperatures. When heated, toxic vapours may be formed.  
 Ingestion: Corrosive. Even small amounts may cause serious damage.

## 12 Ecological information

- **Toxicity** Harmful to aquatic life with long lasting effects.
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** The degradability of the product has not been stated.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Ecotoxicological effects:**
- **Remark:** Harmful to fish
- **Additional ecological information:**
- **General notes:**  
 Water hazard class 2 (German Regulation) (Assessment by list): hazardous for water  
 Do not allow product to reach ground water, water course or sewage system.  
 Must not reach sewage water or drainage ditch undiluted or unneutralized.  
 Danger to drinking water if even small quantities leak into the ground.  
 Harmful to aquatic organisms
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

## 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**  
 Must not be disposed together with household garbage. Do not allow product to reach sewage system.
- **Uncleaned packaging:**
- **Recommendation:**  
 Dispose of waste and residues in accordance with local authority requirements. Waste is classified as hazardous waste. Note that fully cured material is not considered as hazardous waste.  
 Waste from residues: EWC-code: 16 05 08
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

## 14 Transport information

- |                                  |                           |
|----------------------------------|---------------------------|
| · <b>UN-Number</b>               |                           |
| · <b>ADR, IMDG, IATA</b>         | UN2259                    |
| · <b>UN proper shipping name</b> |                           |
| · <b>ADR</b>                     | 2259 TRIETHYLENETETRAMINE |
| · <b>IMDG, IATA</b>              | TRIETHYLENETETRAMINE      |

*(Contd. on page 7)*

# Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

**Trade name: EpoFix Hardener 130ml**

(Contd. of page 6)

- **Transport hazard class(es)**

- **ADR, IMDG, IATA**



- **Class** 8 Corrosive substances.
- **Label** 8

- **Packing group**

- **ADR, IMDG, IATA** II

- **Environmental hazards:**

- **Marine pollutant:** No

- **Special precautions for user**

Warning: Corrosive substances.

- **Danger code (Kemler):**

80

- **EMS Number:**

F-A,S-B

- **Segregation groups**

Alkalis

- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code**

Not applicable.

- **Transport/Additional information:**

- **ADR**

- **Limited quantities (LQ)** 1L
- **Transport category** 2
- **Tunnel restriction code** E

- **UN "Model Regulation":**

UN2259, TRIETHYLENETETRAMINE, 8, II

## 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **National regulations:**

- **Information about limitation of use:**

As a general rule, persons under 18 years of age are not allowed to work with this product. Users must be carefully instructed in the proper work procedure, the dangerous properties of the product and the necessary safety instructions.

- **Other regulations, limitations and prohibitive regulations**

(Contd. on page 8)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 20.08.2014

Revision: 20.08.2014

**Trade name: EpoFix Hardener 130ml***(Contd. of page 7)*

Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive

76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC, including amendments. Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 with amendments. The Chemicals (Hazard Information and Packaging for Supply) Regulations 2009

(S.I 2009 No. 716). The Control of Substances Hazardous to Health Regulations 2002 (S.I 2002 No. 2677) with amendments. The Management of Health and Safety at Work Regulations 1999 (SI 1999 No. 3242), with amendments. The List of Wastes (England) (Amendment) Regulations 2005. (SI 2005 No. 895).

- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge and should assist the user with the safe handling of this material when properly applied. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Department issuing MSDS:** Sales department

- **Contact:**

sales@agarscientific.com

Tel: +44 (0) 1279 813 519

- **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

GB