

# Non-Woven 3-Ply Disposable Looped Face Mask EN14683 Type IIR

(box of 50)



- Our 3-Ply medical Face Mask provides a physical barrier to dust, dirt, debris, smoke, Anti-virus and powder. Multi-layer design blocks air pollution to protect user.
- Folds expand to offer full coverage over nose and mouth.
- Nose piece along top of mask allow for a custom fit.
- Lightweight with elastic earloops, it is comfortable and easy to put on or take off.
- 3ply, non woven PP+meltblown filter+non woven PP
- Size: 17.5\*9.5cm for adults, 14\*9.0cm for children
- Style: with soft flexible ear-loop, fixable nose-clip

## PRODUCT INFORMATION

Ply	3ply (Three layers) 25+25+25
Size	17.5 x 9.5cm
Material	Non-Woven
Certificate	CE & EN14683 Standard
Colour	Blue
Bacterial Filtration	BFE ≥98%
Feature	Convenient, Breathable, Water-proof, Anti-virus

## PACKAGING

Packaging	50pcs per box
-----------	---------------

