

С

ar	Safety data sheet according to 1907/2006/EC, Article	Page 1/6
Date Printed 17.06.2022	Version number 1	Revision Date 17.06.2022
SECTION 1: Ident undertaking	tification of the substance/mixt	ure and of the company/
· 1.1 Product identifier · Trade name: Liquid I	L 2	
No further relevant infor	I uses of the substance or mixture and u	-
• 1.3 Details of the supp • Supplier. Agar Scientific Ltd Parsonage Lane Stansted CM24 8GF United Kingdom sales@agarscientific.co Tel: +44 (0) 1279 813 5		
	btainable from: Technical Support one number: 24 hours: +44 (0)1856 40733	3
SECTION 2: Hazar	ds identification	
Classification accordi	ne substance or mixture ng to Regulation (EC) No 1272/2008 sified, according to the CLP regulation.	
 2.2 Label elements Labelling according to Hazard pictograms Vo 	Regulation (EC) No 1272/2008 Void	

- · Signal word Void
- · Hazard statements Void
- · 2.3 Other hazards
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- vPvB: Not applicable.

SECTION 3: Composition/information on ingredients

- · 3.2 Chemical characterisation: Mixtures
- · Description: Mixture of substances listed below with nonhazardous additions.
- · Dangerous components: Void
- Additional information: For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

- · 4.1 Description of first aid measures
- · General information: No special measures required.
- After inhalation: Supply fresh air; consult doctor in case of complaints.
- After skin contact: Generally the product does not irritate the skin.
- After eye contact: Rinse opened eye for several minutes under running water.
- · After swallowing: If symptoms persist consult doctor.
- 4.2 Most important symptoms and effects, both acute and delayed May cause slight eye irritation (abrasive particles)

(Contd. on page 2)

GB



Safety data sheet according to 1907/2006/EC, Article 31

Version number 1

Revision Date 17.06.2022

Trade name: Liquid L2

Date Printed 17.06.2022

(Contd. of page 1) • **4.3 Indication of any immediate medical attention and special treatment needed** No further relevant information available.

SECTION 5: Firefighting measures

· 5.1 Extinguishing media

· Suitable extinguishing agents:

Use fire extinguishing methods suitable to surrounding conditions.

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide

- · 5.2 Special hazards arising from the substance or mixture Carbon oxides
- 5.3 Advice for firefighters

· Protective equipment: Wear self-contained breathing apparatus for firefighting if necessary

SECTION 6: Accidental release measures

• 6.1 Personal precautions, protective equipment and emergency procedures Use personal protective equipment. Avoid breathing vapours, mist or gas. Ensure adequate ventilation.

- **6.2 Environmental precautions:** Dilute with plenty of water. Do not let product enter drains
- 6.3 Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- 6.4 Reference to other sections
 See Section 7 for information on safe handling.
 See Section 8 for information on personal protection equipment.
 See Section 13 for disposal information.

SECTION 7: Handling and storage

· 7.1 Precautions for safe handling No special measures required.

· Information about fire - and explosion protection: No special measures required.

· 7.2 Conditions for safe storage, including any incompatibilities

· Storage:

· Requirements to be met by storerooms and receptacles:

Avoid temperatures below 2°C. Frost makes the product unusable Avoid storing in steel containers

Store in a cool place. Keep container tightly closed in a dry and well-ventilated place

· Information about storage in one common storage facility: Not required.

- · Further information about storage conditions: None.
- 7.3 Specific end use(s) No further relevant information available.

SECTION 8: Exposure controls/personal protection

- · 8.1 Control parameters
- Additional information about design of technical facilities: No further data; see item 7. • Ingredients with limit values that require monitoring at the workplace:
- The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- Additional information: The lists valid during the making were used as basis.

(Contd. on page 3)

GB ·

Safety data sheet according to 1907/2006/EC, Article 31

Version number 1

Revision Date 17.06.2022

Trade name: Liquid L2

Date Printed 17.06.2022

· 8.2 Exposure controls

- · Personal protective equipment:
- General protective and hygienic measures:
- The usual precautionary measures are to be adhered to when handling chemicals.
- Respiratory protection: Not required.
- · Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Nitrile rubber (thickness 0.4mm) - break through time > 10 min

Eye protection: Goggles recommended during refilling

9.1 Information on basic physical and c General Information	hemical properties
Appearance:	
Form:	Liquid
Colour:	White
Odour:	Characteristic
Odour threshold:	Not determined.
pH-value at 20 °C:	7.5-9.5
Change in condition	
Melting point/freezing point: Initial boiling point and boiling range:	1 °C 100 °C
Flash point:	Not applicable.
Flammability (solid, gas):	Not applicable.
Decomposition temperature:	Not determined.
Auto-ignition temperature:	Product is not selfigniting.
Explosive properties:	Product does not present an explosion hazard.
Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.
Vapour pressure:	Not determined.
Density at 20 °C:	1.3-2.5 g/cm ³

(Contd. of page 2)



Safety data sheet according to 1907/2006/EC, Article 31

Date Printed 17.06.2022

Version number 1

Revision Date 17.06.2022

Trade name: Liquid L2

	(Contd. of page 3)	
· Relative density	Not determined.	
· Vapour density	Not determined.	
· Evaporation rate	Not determined.	
· Solubility in / Miscibility with		
water:	Fully miscible.	
· Partition coefficient: n-octanol/water:	Not determined.	
· Viscosity:		
Dynamic:	Not determined.	
Kinematic:	Not determined.	
· Solvent content:		
VOC (EC)	0.00 %	
· 9.2 Other information	No further relevant information available.	

SECTION 10: Stability and reactivity

- · 10.1 Reactivity No further relevant information available.
- 10.2 Chemical stability
- Thermal decomposition / conditions to be avoided:
- No decomposition if used according to specifications.
- 10.3 Possibility of hazardous reactions No dangerous reactions known.
- · 10.4 Conditions to avoid Avoid frost. Avoid light
- · 10.5 Incompatible materials: Oxidising agents, acids and bases
- 10.6 Hazardous decomposition products: Carbon dioxide and carbon monoxide

SECTION 11: Toxicological information

- · 11.1 Information on toxicological effects
- Acute toxicity Based on available data, the classification criteria are not met.
- · Specific symptoms in biological assay:
- · Skin corrosion/irritation Based on available data, the classification criteria are not met.
- Serious eye damage/irritation Mild eye irritation (abrasive)
- · Respiratory or skin sensitisation Based on available data, the classification criteria are not met.
- · Additional toxicological information:
- · CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)
- Germ cell mutagenicity Based on available data, the classification criteria are not met.
- · Carcinogenicity Based on available data, the classification criteria are not met.
- Reproductive toxicity Based on available data, the classification criteria are not met.
- STOT-single exposure Based on available data, the classification criteria are not met.
- STOT-repeated exposure Based on available data, the classification criteria are not met.
- Aspiration hazard Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- · 12.1 Toxicity
- Aquatic toxicity: No further relevant information available.
- 12.2 Persistence and degradability No further relevant information available.
- · 12.3 Bioaccumulative potential No further relevant information available.
- 12.4 Mobility in soil No further relevant information available.
- Additional ecological information:
- · General notes: Not known to be hazardous to water.

(Contd. on page 5)

agar scientific

Safety data sheet according to 1907/2006/EC, Article 31

Date Printed 17.06.2022

Version number 1

Revision Date 17.06.2022

Trade name: Liquid L2

· 12.5 Results of PBT and vPvB assessment

• **PBT:** Not applicable.

• **vPvB:** Not applicable.

12.6 Other adverse effects No further relevant information available.

SECTION 13: Disposal considerations

· 13.1 Waste treatment methods

· Recommendation Smaller quantities can be disposed of with household waste.

• **Recommendation:** Disposal must be made according to official regulations.

· Recommended cleansing agents: Water, if necessary together with cleansing agents.

SECTION 14: Transport information		
· 14.1 UN-Number · ADR, IMDG, IATA	Void	
 14.2 UN proper shipping name ADR, IMDG, IATA 	Void	
 14.3 Transport hazard class(es) 		
· ADR, ADN, IMDG, IATA · Class	Void	
 14.4 Packing group ADR, IMDG, IATA 	Void	
· 14.5 Environmental hazards:	Not applicable.	
 14.6 Special precautions for user 	Not applicable.	
 14.7 Transport in bulk according to Annex of Marpol and the IBC Code 	II Not applicable.	
· UN "Model Regulation":	Void	

SECTION 15: Regulatory information

 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I None of the ingredients is listed.
- 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge and should assist the user with the safe handling of this material when properly applied. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- · Department issuing SDS: Sales department
- · Contact:

sales@agarscientific.com Tel: +44 (0) 1279 813 519

• Abbreviations and acronyms: ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

(Contd. on page 6)

GB –



Safety data sheet according to 1907/2006/EC, Article 31

Date Printed 17.06.2022

Version number 1

Revision Date 17.06.2022

Trade name: Liquid L2

	(Contd. of page 5
IMDG: International Maritime Code for Dangerous Goods	
IATA: International Air Transport Association	
GHS: Globally Harmonised System of Classification and Labelling of Chemicals	
EINECS: European Inventory of Existing Commercial Chemical Substances	
ELINCS: European List of Notified Chemical Substances	
CAS: Chemical Abstracts Service (division of the American Chemical Society)	
VOC: Volatile Organic Compounds (USA, EU)	
PBT: Persistent, Bioaccumulative and Toxic	
vPvB: very Persistent and very Bioaccumulative	
* Data compared to the previous version altered.	