

Ferritin 6X Crystallised

AGR1200

(CADMIUM REMOVED)

Source:	Horse Spleen
Total Solids:	50.0 mg/ml
Ferritin Solids:	49.8 mg/ml
Absorption: A440nm	44. mg/ml
Protein: By Nitrogen Analyser:	5.4 g%
pH (Direct):	7.6
Iron:	8.7 mg/ml
Cadmium:	40.0 ug/ml
Sterility:	Bacteriologically sterile by U.S.P X1X procedures.
Electrophoresis:	Cellulose acetate: Barbital Buffer, pH 9.6, 0.075 ionic strength
Purity:	95%