

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name: CitiFluor CFM-3 Mountant Solution**
- **Article number:** AGR1332-10
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
A buffered aqueous-glycerol solution of a phenolic antioxidant for use in fluorescence microscopy
- **Application of the substance / the preparation:** No further relevant information.
- **1.3 Details of the supplier of the safety data sheet**
- **Supplier.**
Agar Scientific Ltd
Parsonage Lane
Stansted CM24 8GF
United Kingdom
sales@agarscientific.com
Tel: +44 (0) 1279 813 519
- **Further information obtainable from:** Technical Support
- **1.4 Emergency telephone number:** 24 hours: +44 (0)1856 407333

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**



Acute Tox. 4 H332 Harmful if inhaled.
 Eye Irrit. 2 H319 Causes serious eye irritation.
 Skin Sens. 1 H317 May cause an allergic skin reaction.

- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**
The product is classified and labelled according to the CLP regulation.
- **Hazard pictograms**



GHS07

- **Signal word** Warning
- **Hazard-determining components of labelling:**
pyridine 1-oxide
propyl gallate
- **Hazard statements**
H332 Harmful if inhaled.
H319 Causes serious eye irritation.
H317 May cause an allergic skin reaction.
- **Precautionary statements**
P101 If medical advice is needed, have product container or label at hand.
P102 Keep out of reach of children.
P103 Read carefully and follow all instructions.
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves / eye protection / face protection.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

Trade name: CitiFluor CFM-3 Mountant Solution
(Contd. of page 1)

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P312 Call a POISON CENTER/doctor if you feel unwell.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

· 3.2 Chemical characterisation: Mixtures

- **Description:** Mixture of substances listed below with nonhazardous additions.

· **Dangerous components:**

CAS: 694-59-7	pyridine 1-oxide	50.0%
EINECS: 211-774-6	⚠ Acute Tox. 4, H312; Acute Tox. 4, H332; Eye Irrit. 2, H319	
CAS: 56-81-5	glycerol	40.0%
EINECS: 200-289-5	substance with a Community workplace exposure limit	
CAS: 121-79-9	propyl gallate	5.0%
EINECS: 204-498-2	⚠ Acute Tox. 4, H302; Skin Sens. 1, H317	

- **Additional information:** For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

· 4.1 Description of first aid measures

· **General information:**

Immediately remove any clothing soiled by the product.

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

· **After inhalation:**

Supply fresh air.

No adverse effects are anticipated from inhalation.

Move the exposed person to fresh air at once. Rinse nose and mouth with water. Get medical attention if any discomfort continues. Rinse and clean nose and mouth with water. Artificial respiration may be administered by suitably qualified first aiders.

· **After skin contact:**

Wash with water and soap and rinse thoroughly.

Remove affected person from source of contamination. Remove contaminated clothing. Wash the skin immediately with soap and running water. Take special care to clean folds, crevices and groin. Get medical attention.

· **After eye contact:**

Rinse opened eye under running water. If symptoms persist, consult a doctor.

Check for contact lenses which must be removed from the eyes before rinsing. Promptly rinse eyes with plenty of clean water while lifting the eyelids. Continue to rinse for at least 15 minutes. Continue until the eyes are free of all traces of contamination. Get medical attention if any discomfort or irritation persists.

· **After swallowing:**

NEVER MAKE AN UNCONSCIOUS PERSON VOMIT OR DRINK FLUIDS! Get medical attention if any discomfort continues. If patient vomits keep head low to prevent vomit entering the lungs.

· **4.2 Most important symptoms and effects, both acute and delayed**

No further relevant information available.

(Contd. on page 3)

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

Trade name: CitiFluor CFM-3 Mountant Solution

(Contd. of page 2)

- **4.3 Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
- **Suitable extinguishing agents:**
Use fire extinguishing methods suitable to surrounding conditions.
This product is not flammable. Use fire extinguishing media appropriate for surrounding materials
- **For safety reasons unsuitable extinguishing agents:**
Do not use water jet as an extinguisher as this may spread the fire.
- **5.2 Special hazards arising from the substance or mixture**
During fire, toxic gases (CO, CO₂, NO_x) may be formed.
- **5.3 Advice for firefighters**
- **Protective equipment:**
Positive pressure self-contained breathing apparatus.
Wear self contained breathing apparatus and full protective clothing. Keep all unnecessary people away. fire water run-off must not be allowed to contaminate ground or enter drains, sewers or water courses. Provide bunding against fire water run-off.

SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**
Wear protective gloves. Eye contact MUST be prevented by means of suitable personal protection equipment
- **6.2 Environmental precautions:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
Do not allow ANY environmental contamination. Prevent discharge into drains, water courses or onto the ground.
- **6.3 Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Dispose contaminated material as waste according to item 13.
Ensure adequate ventilation.
Stop leak if possible without risk. Absorb in vermiculite, dry sand or earth and place into containers.
Flush with plenty of water to clean spillage area. Flushings must not be allowed to contaminate ground or enter drains, sewers or water courses.
- **6.4 Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
Avoid spilling, skin and eye contact. Wash thoroughly after handling.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
Store in tightly closed container in a dry, cool and well-ventilated place. Keep away from food, drink and animal feeding stuffs.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.

(Contd. on page 4)

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

Trade name: CitiFluor CFM-3 Mountant Solution

(Contd. of page 3)

- **Further information about storage conditions:** Keep container tightly sealed.
- **7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
- **Additional information about design of technical facilities:**
Provide adequate ventilation. Observe occupational exposure limits and minimise the risk of inhalation of vapours.

Ingredients with limit values that require monitoring at the workplace:

56-81-5 glycerol	
WEL	Long-term value: 10 mg/m ³

- **Additional information:** The lists valid during the making were used as basis.
- **8.2 Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
Keep away from foodstuffs, beverages and feed.
Remove all soiled and contaminated clothing
Wash hands before breaks and at the end of work.
Avoid contact with the eyes.
Avoid contact with the eyes and skin.
DO NOT SMOKE IN WORK AREA! Wash hands at the end of each work shift and before eating, smoking and using the toilet. Wash promptly if skin becomes wet or contaminated. Promptly remove any clothing that becomes contaminated. Use appropriate skin cream to prevent drying of skin. When using do not eat, drink or smoke.
- **Respiratory protection:**
In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.
In case of inadequate ventilation of when the product is heated use suitable respiratory equipment with gas filter (type A2).
- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Wear suitable protective gloves conforming to EN 374. Seek recommendation from manufacturer or supplier. After using gloves the hands should be washed and dried thoroughly and a suitable moisturiser applied.

- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 5)

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

Trade name: CitiFluor CFM-3 Mountant Solution

(Contd. of page 4)

· Eye protection:



Tightly sealed goggles

Use approved safety goggles or face shield.

 · **Body protection:** Wear appropriate clothing to prevent any possibility of skin contact.

SECTION 9: Physical and chemical properties

 · **9.1 Information on basic physical and chemical properties**

 · **General Information**

 · **Appearance:**

Form:	Liquid
Colour:	Colourless
· Odour:	Odourless
· Odour threshold:	Not determined.

· pH-value: Not determined.

 · **Change in condition**

Melting point/freezing point:	Undetermined.
Initial boiling point and boiling range:	Undetermined.

· Flash point: Not applicable.

· Flammability (solid, gas): Not applicable.

· Decomposition temperature: Not determined.

· Auto-ignition temperature: Product is not selfigniting.

· Explosive properties: Product does not present an explosion hazard.

 · **Explosion limits:**

Lower:	Not determined.
Upper:	Not determined.

· Vapour pressure: Not determined.

· Density: Not determined.

· Relative density: Not determined.

· Vapour density: Not determined.

· Evaporation rate: Not determined.

 · **Solubility in / Miscibility with water:**

Fully miscible.

· Partition coefficient: n-octanol/water: Not determined.

 · **Viscosity:**

Dynamic:	Not determined.
Kinematic:	Not determined.

 · **Solvent content:**

Organic solvents:	40.0 %
VOC (EC)	0.00 %

Solids content: 0.0 %

 · **9.2 Other information**

Refractive index 1.52 @20°C

GB

(Contd. on page 6)

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

Trade name: **CitiFluor CFM-3 Mountant Solution**

(Contd. of page 5)

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability** Stable under normal conditions of storage and use.
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** Strong oxidising agents
- **10.6 Hazardous decomposition products:**
Decomposition can lead to the formation of toxic gases or fumes, including carbon monoxide (CO), carbon dioxide (CO₂) and nitrous gases (NO_x).

SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
- **Acute toxicity**
Harmful if inhaled.
- **Specific symptoms in biological assay:**
- **Skin corrosion/irritation**
May cause defatting of the skin, but is not an irritant. Not a skin sensitiser.
- **Serious eye damage/irritation**
Causes serious eye irritation.
Irritating and may cause redness and pain
- **Respiratory or skin sensitisation**
May cause an allergic skin reaction.
- **Other information (about experimental toxicology):**
Inhalation: Not relevant at room temperatures. When heated, irritating vapours may be formed.
Ingestion: May cause discomfort if swallowed.
- **Additional toxicological information:**
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- **12.1 Toxicity**
The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** This product is expected to be biodegradable
- **12.3 Bioaccumulative potential**
This product does not contain any substances expected to be bioaccumulating
- **12.4 Mobility in soil** The product is completely water-soluble and may spread in water systems.
- **Other information:**
The product may contribute to an excessive enrichment of the aquatic environment with nutrients.
- **Additional ecological information:**
- **General notes:**
Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water
Do not allow product to reach ground water, water course or sewage system.
Danger to drinking water if even small quantities leak into the ground.

(Contd. on page 7)

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

Trade name: CitiFluor CFM-3 Mountant Solution
(Contd. of page 6)

- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Other adverse effects** No further relevant information available.

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**
Must not be disposed together with household garbage. Do not allow product to reach sewage system.
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

SECTION 14: Transport information

- | | |
|---|-----------------|
| <ul style="list-style-type: none"> · 14.1 UN-Number · ADR, ADN, IMDG, IATA | Void |
| <ul style="list-style-type: none"> · 14.2 UN proper shipping name · ADR, ADN, IMDG, IATA | Void |
| <ul style="list-style-type: none"> · 14.3 Transport hazard class(es) · ADR, ADN, IMDG, IATA · Class | Void |
| <ul style="list-style-type: none"> · 14.4 Packing group · ADR, IMDG, IATA | Void |
| <ul style="list-style-type: none"> · 14.5 Environmental hazards: · Marine pollutant: | No |
| <ul style="list-style-type: none"> · 14.6 Special precautions for user | Not applicable. |
| <ul style="list-style-type: none"> · 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code | Not applicable. |
| <ul style="list-style-type: none"> · UN "Model Regulation": | Void |

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge and should assist the user with the safe handling of this material when properly applied. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**
- H302 Harmful if swallowed.
- H312 Harmful in contact with skin.
- H317 May cause an allergic skin reaction.

(Contd. on page 8)

Safety data sheet

according to 1907/2006/EC, Article 31

Date Printed 17.02.2022

Version number 1

Revision Date 17.02.2022

Trade name: CitiFluor CFM-3 Mountant Solution*(Contd. of page 7)*

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

· **Department issuing SDS:** Sales department· **Contact:**

sales@agarscientific.com

Tel: +44 (0) 1279 813 519

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity – Category 4

Eye Irrit. 2: Serious eye damage/eye irritation – Category 2

Skin Sens. 1: Skin sensitisation – Category 1

· *** Data compared to the previous version altered.**

GB