

X-Ray Standards for Transmission Electron Microscopy


AGS1981

These standards are supplied as fine powders dispersed onto holey carbon films and are selected from a range of certified natural and synthetic compounds.

Normally the holey carbon films are supported by 400 mesh copper grids although alternative grid materials can be supplied.

The rare earth set contains 14 compound standards:

LaF₃ SmF₃ DyF₃ YbF₃
CeAl₂ EuF₃ HoF₃ LuF₃
PrF₃ GdF₃ ErF₃
NdF₃ TbF₃ TmF₃



AGS1981 Rare earth TEM X-ray standard set.