

X-Ray Standards for Transmission Electron Microscopy

AGS1981

These standards are supplied as fine powders dispersed onto holey carbon films and are selected from a range of certified natural and synthetic compounds.

Normally the holey carbon films are supported by 400 mesh copper grids although alternative grid materials can be supplied.

The rare earth set contains 14 compound standards:

LaF_3	SmF_3	DyF_3	YbF_3
CeAl_2	EuF_3	HoF_3	LuF_3
PrF_3	GdF_3	ErF_3	
NdF_3	TbF_3	TmF_3	

S1981 Rare earth TEM X-ray standard set