


X-Ray Standards for Transmission Electron Microscopy

AGS1982

Metal foil standards. Set of 25 AGS1982

For STEM applications there is a thin foil standard set consisting of 25 metal foils each measuring 3mm diameter and 0.1mm thick.

Aluminium	Magnesium	Tantalum
Cadmium	Molybdenum	Tin
Cobalt	Nickel	Titanium
Copper	Niobium	Tungsten
Gold	Palladium	Vanadium
Hafnium	Platinum	Yttrium
Indium	Rhodium	Zinc
Iron	Silver	Zirconium
AISI 316 S/S		



Alternatively, these foils can be supplied as single standards minimum order 5: See AGS1983