

X-Ray Standards for Transmission Electron Microscopy

AGS1983

Metal foil standards. Set of 5 AGS1983

For STEM applications there is a thin foil standard set consisting of 5 metal foils each measuring 3mm diameter and 0.1mm thick.

Please choose from list below:

Aluminium	Magnesium	Tantalum
Cadmium	Molybdenum	Tin
Cobalt	Nickel	Titanium
Copper	Niobium	Tungsten
Gold	Palladium	Vanadium
Hafnium	Platinum	Yttrium
Indium	Rhodium	Zinc
Iron	Silver	Zirconium
AISI 316 S/S		

